


(GB)

BODY COMPOSITION ANALYZER SC-330

Instruction manual



(GB)

Table of Contents

Before use
Safety Notes Part names / connecting method Various settings 10
How to use
Operating Instructions 24 When using as a body composition monitor 24 Target body fat ratio 34 When using as a scale 35 Various criteria 36
When necessary
Troubleshooting 42 Connection with a personal computer 44 Technical Notes 48 New regression equation of BMR
(basal metabolic rate)50 Specifications54
Specifications

Applications

- This equipment can be used in the screening of certain adult diseases and conditions related body weight and composition.
- It can be used in the monitoring and prevention of conditions caused by excessive deposits of fat tissue, such as diabetes, hyperlipidemia, cholelithiasis and fatty liver.
- It can be used in the monitoring of changes in individuals body composition, related to differences in the ratio of fat tissue to lean.
- It can be used to assess the effectiveness of individuals nutrition and exercise programmers, both for health and physical fitness.

Benefits

- 1. This product is simple to use, and requires no specialized facilities or expertise to take measurements.
- 2. Measurements can be taken quickly and easily, causing minimal inconvenience to the patient during measurement.

Before use (Safety Notes)

Safety Notes

Caution Symbols

Thank you for purchasing this precision crafted Tanita product. For optimum performance and safety, please familiarize yourself with the Caution Symbols below. These symbols are designed to alert the user to potential hazards when using this equipment. Ignoring these Caution Symbols may result in serious injury, or damage to the product.

Please be sure to review before proceeding with the INSTRUCTION MANUAL.

WARNING This symbol indicates the possibility of serious injury if the product is mishandled or instructions are ignored.



CAUTION This symbol indicates the possibility of physical injury or equipment damage if instructions are ignored.



This symbol indicates general precautions that should be taken when using this product.



Individuals with a Pacemaker or Other Internal Medical Devices

This equipment sends a weak electrical current through the body during measurement. Individuals who have internally implanted medical devices, such as Pacemakers, should not use this equipment due to the risk of malfunction to the device that may be caused by the weak electrical current.

Inserting and Removing the Power Cord

To reduce the risk of electric shock or product damage, never insert or remove the power cord with

- Do not under any circumstances dismantle or alter the device, as this could result in electric shock or injury as well as adversely affect the precision of measurements.
- To prevent fire hazard

Use only a correctly wired (100-240VAC) outlet, and do not use a multiple outlet extension cable.

Measurements for physically disabled persons

Physically disabled persons should not attempt to take measurements alone, but instead should have their caretakers assist them in using the device.



Cross Contamination

The Body Composition Analyzer should be used with bare feet. Please be sure to clean the scale platform with appropriate disinfectant after each use. Never pour any liquid directly on the scale platform, as it may leak and cause internal damage. Use a soft cloth and appropriate ethyl alcohol to wipe off platform. Do not wipe the platform with strong chemicals.

The data provided by this machine, as well as any supplementary information such as diet or exercise programs based on this data, should be interpreted by a licensed professional.

Please make sure you place the Weighing Platform on a level and stable surface. If the equipment is used when the Weighing Platform is unstable because not all feet are on the surface, there may be a risk of stumbling or inaccurate measurement.

Never jump on the Weighing Platform, there may be a risk of stumbling and malfunction of the equipment.

When handling printer unit, avoid any sharp edges.

For the SC-330; Ensure you use the original AC adapter (MODEL: SA165A-0950U-3). Using an AC adapter other than the original one may cause malfunction.

Do not insert or remove the plug by the cable.

Maintenance

Since this equipment is accurately manufactured and adjusted, please observe the following instructions.

- Never disassemble the equipment as this may cause malfunction. Users must not disassemble or adjust this equipment. Please inspect the equipment in accordance with the regulations in your country.
- Unplug the unit from the wall outlet when it will not be in use for long periods of time.
- In order to reduce the risk of a short circuit, please keep any liquid or metal objects (paper clips, etc.) away from the printer.
- Keep the electrodes clean by wiping them with disinfectant.
- Do not drop the unit, and avoid locations with constant vibration.
- Do not put this equipment in direct sunlight, close to heaters or near direct draughts from air conditioners.
- When transferred to any location where there is a difference of more than 20 degrees centigrade (40 degrees Fahrenheit), wait 2 hours before using.
- When disposing of this unit, please do so in accordance with the prevailing regulations in each country.

General Instructions for Accurate Measurement

This equipment sends out a very weak electric current to measure impedance (electrical resistance) of the body. Therefore, in principle, users need to use this equipment with bare feet. Moreover, since impedance fluctuates in accordance with the distribution of body fluid, please observe the following instructions for accurate measurement.

- To prevent a possible discrepancy in measured values, avoid taking measurements after vigorous exercise until sufficiently rested.
- To prevent inaccurately low body fat percentage measurements and other measurement errors, always hold both arms straight down when taking measurements.
- · As changes in body water and body temperature can have a major impact on measurements, measurements should be made every day at the same time under similar conditions (always urinating before taking measurements, etc.) to get a more accurate picture of the measurements over time.
- · Ensure that your arms are not touching your side and that the inner thighs are not touching each other during measurements; if necessary, place a dry towel between your arm and side and/or between your thighs.
- · Also, make sure the soles of feet are free of excess dirt, as this may also act as a barrier to the mild current.
- False results may be reported after excessive food /fluid intake, or after periods of intense exercise.
- * For further details, see the Technical Notes on page 48.
- This equipment is designed for the majority of the population leading healthy lives with a regular lifestyle. For people suffering from sickness, or whose lifestyle is very different from the norm, it is recommended that the data from this product should not be used as an absolute value, but rather as a reference to observe the rate change.
- * For further details, see the Technical Notes on page 48.
- · Measurement is sometimes impossible on a surface that is strongly vibrating. In this case, please move the equipment onto a surface with little vibration.
- Do not take measurements while using transmitters, such as mobile phones, which may affect readings.

<Usage Conditions>

Temperature Range for Use : 0°C — 35°C

: 30% — 80% (without condensation) **Relative Humidity**

<Storage Conditions>

Temperature Range of Environment :-10°C — 50°C

Range of Relative Humidity : 10% — 90% (without condensation)

To avoid malfunctions, avoid storing the equipment where there is direct sunlight, significant temperature changes, the risk of dampness, a large amount of dust, in the vicinity of fires, or where there is the risk of receiving vibrations or shocks.

<Power Source>

Model Name	SC-330
Frequency Range	50 / 60Hz
Electric Current Range	1.5A


Before use (Product Assembly and Components)

Product Assembly and Components


Remote display version


Control Box Printer cover Operating key

the product for details.



Column mounted version







Make sure you place the Weighing Platform on a stable, level surface. If the Weighing Caution Platform is not stable because not all the feet are on the surface, for example, there is a risk of stumbling or inaccurate measurement.


(4 positions)


Platform

Electrode

(GB)

Before use (Product Assembly and Components)





radio waves, such as mobile phones.

⚠ Warning

controller

•To avoid electric shock, do not insert or remove the plug with wet



- To avoid electric shock, do not use the equipment near water. • To avoid measurement error, do not measure while using equipment that generates
- Use only the original AC adapter (MODEL: SA165-0950U-3). Using AC adapters other than the original one may cause malfunction, smoke or fire.

Symbols and their Meanings

•	Power On	Ċ	Power Off		Direct current	\Longrightarrow	Input, Output		Class II Equipment
FEED	Advances the paper	\triangle	Caution Refer to the attached notes.	Male 🟺	Male	Female 🖡	Female	РТ	Clothes Weight Setting

- Please change the paper roll when red lines appear along the sides of the paper.
- Please pay careful attention to avoid injury from the sharp edge.
 - Please turn off the machine before clearing Paper jams.

Setting of the printer paper roll


Press **O/O** to turn on the power.

• After all lamps light up, the model number is displayed as [330], and [[]]_{kg} is displayed.



(1) Press OPEN.

(2) Remove the printer cover.




(3) Set up the print roll paper

• Remove the adhesive of the printer paper and draw it out approx. 10 cm.

Return the printer cover. • In the case that [[]PEn] is displayed,


⇒The printer cover is open, so please close it again properly (property (property 43).



Press the FEED , and cut off the excessive paper.



• In the case that the automatic cutting is set to "OFF," the automatic cutting does not function (repage 14).



Setting completed.

(GB)

Before use (Preparation)

Call up the setting item selection screen.

Press \bigcirc/\bigcirc to turn on the power.

• After all lamps light up, the model number is displayed, and [[]] sq. is displayed.



Press Set Up

(GB)

Before use (Various settings)

• The setting item input screen is displayed.

Note • 39 – 67 are the setting items related to print out items (page 20).

• When the various settings are all completed,

⇒Press Set Up on the "Setting item selection" screen (it returns to the tare input screen).

Setting items

1	Date and time(page 12)
2	Number of sheets to print
	Body composition mode (F page 13)
3	Number of sheets to print
	Weight only mode (🎓 page 13)
4	Auto Cutting the printer paper (page 14)
5	Beep sound (page 14)
6	Display Fat % Healthy range (page 15)
7	ID No. (page 15)
8	Measurement flow (page 16)
9	Athletic mode (page 16)
10	Input unit of height (page 17)
11	Automatic determination time (page 17)
18	Target body fat ratio (page 18)
19	Select language (FF page 18)
20	Print item preset (page 18)

Print item setting ON / OFF

39	TANITA Logo (FF page 22)
40	Category name (page 22)
41	Date (page 22)
42	Serial number (page 22)
43	Memo space (page 22)
44	ID No (page 22)
45	Fat mass (Fat page 22)
46	Fat free mass (FF page 22)
47	Muscle mass (page 22)
48	Total body water (🎓 page 22)
49	Total body water % (🖙 page 22)
50	Bone mass (page 22)
51	BMR (ref page 22)
52	Metabolic age (page 22)
53	Visceral fat rating (page 22)
54	BMI (F page 22)
55	The Rohrer's index (page 22)
56	Ideal body weight (F page 22)
57	Degree of obesity (page 22)
58	Desirable range Body fat (page 22)
59	Graph Fat % (page 22)
60	Graph BMI (F page 22)
61	Graph Visceral fat rating (page 22)
62	Graph Muscle mass (FF page 22)
63	Graph BMR (F page 22)
64	Physique rating (r page 22)

Note


- *When various settings are continuously carried out, ⇒press each number to set.
- *The set contents are memorized until they are changed next time.


⇒press Set Up on the "setting item selection screen" (it returns to the tare input screen).

^{*}When various settings are all completed,

GB)

Before use
(Various settings)





Note

- *When various settings are continuously carried out,

 ⇒press each number to set.
- *The set contents are memorized until they are changed next time.

press Enter / Next

press Set Up on the "setting item selection screen" (it returns to the tare input screen).

After inputting the numeric value,

• It returns to the "setting item selection" screen.

Body Type

Set ON / OFF of the automatic cutting of the printer paper (Continued from page 10).

and then press Enter / Next

• The "ON or OFF of automatic cutting" selection screen is displayed.



Before use (Various settings)

Set ON or OFF of the automatic cutting. • The default is "0.off." ("1. on" for valid, "0. off" for

> • To correct the input value, or cancel the input, ⇒press **CE** (the input is deleted).



After inputting the numeric value, press (Enter / Next).

• It returns to the "setting item selection" screen.



Set ON / OFF of the beep sound (Continued from page 10).

and then press Enter / Next).



Set ON or OFF of the beep sound.

⇒press **CE** (the input is deleted).

• The default is "1.on." ("1. on" for valid, "0. off" for • To correct the input value, or cancel the input,



After inputting the numeric value, press (Enter / Next)

• It returns to the "setting item selection" screen.



Set ON / OFF of the Fat % Healthy range display (Continued from page 10).

Press 6 and then press Enter / Next

• The "ON or OFF of the Fat % Healthy range" selection screen is displayed.



Set ON or OFF of the Fat % Healthy range display when using the body composition



• The default is "1.on." ("1. on" for valid, "0. off" for

• To correct the input value, or cancel the input, \Rightarrow press **CE** (the input is deleted).



After inputting the numeric value, press Enter / Next

• It returns to the "setting item selection" screen.



Press 7 and then press Enter / Next

• The "with or without an ID" setting screen is displayed.



Set with or without an ID.

• The default is "0.off." ("1. on" for valid, "0. off" for



• To correct the input value, or cancel the input, \Rightarrow press **CE** (the input is deleted).

After inputting the numeric value, press Enter / Next

• It returns to the "setting item selection" screen.



Note

- *When various settings are continuously carried out, press each number to set.
- *The set contents are memorized until they are changed next time.

- *When various settings are all completed,
- ⇒press Set Up on the "setting item selection screen" (it returns to the tare input screen).

Select the measurement flow (Continued from page 10).

Press 8 and then press Enter / Next

• The "ON or OFF of the one step mode" selection screen is displayed.

• The one step mode is

⇒a mode to measure body weight after inputting personal information.

Body Type

Set ON or OFF of the one step mode when using the body composition analyzer.



Before use (Various settings)

- The default is "0. off." ("1. ON" for valid, "0. off" for invalid.)
- To correct the input value, or cancel the input, \Rightarrow press **CE** (the input is deleted).

After inputting the numeric value, press Enter / Next

• It returns to the "setting item selection" screen.



Set ON or OFF of the athlete selection (Continued from page 10).

Press 9 and then press Enter / Next .

• The "Set ON / OFF of the athletic mode selection" setting screen is displayed.



Set ON or OFF of the athletic mode selection when using the body composition analyzer.

- The default is "1. on." ("1. on" for valid, "0. off" for
- To correct the input value, or cancel the input, \Rightarrow press (CE) (the input is deleted).



• It returns to the "setting item selection" screen.



Elén

Note

- *When various settings are continuously carried out, press each number to set.
- *The set contents are memorized until they are changed next time.

Set the input unit of height (Continued from page 10).

and then press Press 1 Enter / Next

• The "input unit of height" selection screen is displayed.



Selects the input unit of height. • The default is "1. on."

0. off: sets 0.1 cm unit input 1. on: sets 1 cm unit input

• To correct the input value, or cancel the input,

press (CE) (the input is deleted).



After inputting the numeric value, press (Enter / Next)

• It returns to the "setting item selection" screen.

Set the automatic determination time when inputting (Continued from page 10).

and then press Press 1 **Enter / Next**

> • The "automatic determination time when inputting" setting screen is displayed.



- Automatic determination when inputting is,
- ⇒a function to determine the input value automatically even without pressing **Enter / Next**), after inputting the numeric value.

Set the automatic determination time when inputting.



- The default is 5 seconds ("5"). (input range: 0 9). *If "0" is set, it is not automatically determined.
- To correct the input value, or cancel the input, ⇒press (CE) (the input is deleted).





- After inputting the numeric value, press Enter / Next
 - It returns to the "setting item selection" screen.

*When various settings are all completed,

⇒press Set Up on the "setting item selection screen" (it returns to the tare input screen).

Before use (Various settings)





Set the print item preset (Continued from page 10). and then press Press (

Enter / Next

• The "print out preset" setting screen is displayed.



Set the print item preset. • Set with 1 - 3

• The print items preset are \Rightarrow a function to set the print items of the pattern 1 – 3 to ON automatically at once (page 19). The default is "1" (pattern 1).

• When ON or OFF of print items are changed with the "setting of the print out item" after the print out preset (repage 22), the state set for the last time is valid.

After inputting the numeric value, press Enter / Next .

• It returns to the "setting item selection" screen.



[Lists of contents of the print item preset]

D				Body cor	nposition	monitor				
Print item	1	Patte	rn 1	2	Patte	rn 2	3	Patte:	rn 3	Scale
Body type	Standard	Athletic	Child	Standard	Athletic	Child	Standard	Athletic	Child	
TANITA Logo	$\overline{}$	\vee	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	\vee	$\overline{}$	\vee
Category name	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	\vee	$\overline{}$	$\overline{}$	\vee	$\overline{}$	\vee
Date	$\overline{}$	\vee	$\overline{}$	\vee	\vee	$\overline{}$	\vee	\vee	$\overline{}$	\vee
Serial number	$\overline{}$	\vee	$\overline{}$	$\overline{}$						
Memo Space	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	\vee	$\overline{}$	$\overline{}$	\vee	$\overline{}$	$\overline{}$
ID No	$\overline{}$	\vee	$\overline{}$	\vdash	\vee	$\overline{}$	\vee	\vee	$\overline{}$	\vee
Fat mass	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$				
Fat free mass	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$				
Muscle mass	$\overline{}$	\vee	$\overline{}$	\vee	\vee	$\overline{}$				
Total body water	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$				
Total body water %	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$				
Bone mass	$\overline{}$	\vee		\vdash	\vee					
BMR	$\overline{}$	$\overline{}$		\vee	\vee					
Metabolic age	$\overline{}$	\vee		\vee	\vee					
Visceral fat rating	\vee	\vee		V	\vee					
BMI	$\overline{}$	\vee	$\overline{}$	\vee	\vee	$\overline{}$	$\overline{}$	\vee	$\overline{}$	
Rohrer's index										
Ideal body weight	$\overline{}$			V						
Degree of obesity	$\overline{}$			\vdash						
Desirable range	$\overline{}$	\vee	$\overline{}$							
Graph Fat %	$\overline{}$	\vee	$\overline{}$	$\overline{}$	$\overline{}$	$\overline{}$				
Graph BMI	\vee	\vee		\vdash	\vee					
Graph Visceral fat rating	$\overline{}$	$\overline{}$								
Graph Muscle mass	$\overline{}$	$\overline{}$								
Graph BMR	\vee	\vee								
Physique rating										

- The items marked with "

 " are printed.
- The items marked with "

 " can be selected to print (

 page 22)
- * See page 20 for an example of preset print.
- *When various settings are all completed,

⇒press Set Up on the "setting item selection screen" (it returns to the tare input screen).

Note

- * When various settings are continuously carried out, press each number to set.
- *The set contents are memorized until they are changed next time.

-LOGO -Serial No.

Fat %

FFM

TBW

Bone mass*

entire bone.

is higher.

BMI

-Metabolic age*

• The default is 00000001.

Adds 1 each time it measures.

• When it is set with an ID, it is printed

• Fat % is amount of body fat as a

• Fat Free Mass is comprised of muscle,

bone, tissue, water, and all other fat

•Total Body Water is the amount of

water retained in the body. TBW is said

to comprise between 50% - 70% of total

body weight. Generally, men tend to

have higher water weight than women

• Bone mineral amount included in the

• Metabolic age is evaluated young when

a muscular amount is larger, and BMR

• Calculated with "weight (kg) / height

* The standard value is for the Standard mode. In the case of the Athletic mode, the standard value is just a reference. And for those who are 17

years old or younger, only the body fat % is displayed as the standard value. The muscle mass, total body water and the estimated bone mass for those who are 17 years old or

• Desirable Range 18.5 - 24.9

younger are for reference.


due to a greater amount of muscle.

out. (The default is without an ID.)


proportion of body weight.

free mass in the body.


[In the case to select the print item preset "1"]



[In the case to select the print item preset "2"]



[In the case to select the print item preset "3"]



Please consult your doctor before you start a body weight management program. Tanita is not responsible for the target body fat ratio.

Please consult your doctor before you start a body weight management program. Tanita is not responsible for the target body fat ratio.

*18 - 99 yaers only

Before use (Various settings)

(GB)

Set the items to print out (Continued from page 10).

Select the number to set with the numeric keys (3 and press (Enter / Next).

<In the case of selecting 39>

• The setting screen is displayed.

Note

- See the next page for the setting number of each item.
- When the print item preset is set, after setting the items to print out (page 18), the items to print may be changed automatically. Please confirm, "Print items preset list" (repage 19).

Print item setting ON / OFF

39	TANITA Logo
40	Category name
41	Date
42	Serial number
43	Memo space
44	ID No
45	Fat mass
46	Fat free mass
47	Muscle mass
48	Total body water
49	Total body water %
50	Bone mass
51	Basal metabolic rate (BMR)

52	Metabolic age
53	Visceral fat rating
54	BMI
55	Rohrer's index
56	Ideal body weight
57	Degree of obesity
58	Desirable range Body fat %
59	Graph Fat %
60	Graph BMI
61	Graph Visceral fat rating
62	Graph Muscle mass
63	Graph BMR
64	Physique rating

Selects ON or OFF of the item to print < out.

Note

- "0. off" is not to print out and "1. on" is to print out.
- To correct the input value, or cancel the input, \Rightarrow press (CE) (the input is deleted).

After inputting the numeric value, press Enter / Next .

• It returns to the "setting item selection" screen.



Note


- *When various settings are continuously carried out, press each number to set.
- *The set contents are memorized until they are changed next time.

*When various settings are all completed,

⇒press Set Up on the "setting item selection screen" (it returns to the tare input screen).

 \triangle

This explains the procedure when the printer is turned on. Please be aware that the display may be different if the number of print outs is set to ①.



- When the equipment has been transferred to any location where there is a temperature difference of 20°C or more, wait for at least two hours before using it.
- In taking measurements, please keep the person away from the unit, who uses transmitters such as a mobile phone avoid causing margin errors.

About Athletic Mode

- It is recommended for those who are 18 years old or older and meet the following conditions to select "Athletic Mode" and measure as reference values.
- Those who exercise for 12 hours or more per a week.
- Those who belong to a sport team or a sport organization with the aim of participation in competition, etc.
- Those who exercise to build up like a bodybuilder.
- Those who are professional athletes.

Attention

- The posture when measuring
- Stand with both feet parallel on the electrodes.
- Stand fullface without bending knees.
- The age input range is 5 99 years old. Input age 99 for those who are 100 years old or older.

Note

- False results may be reported after excessive food/fluid intake, or after periods of intense exercise.
- When the user is 18 years or older, Athletic mode can be used.
- When Clothes weight is input, Clothes weight is subtracted from measurements and it displays it as Weight.

In the case of the standard flow (not the one step mode)


In the standard flow, after measuring weight, personal data is input and then body composition is measured.

Weight measurement

Personal data input

Body composition measurement

Press **O/O** to turn on the power.



Check that the body composition monitor is selected and input the clothes weight.


Input it by pressing 0 - 9 and .

Note

- The clothes weight (preset tare) can be input in the range of 0.0 10.0 kg.
- To correct the input value,

ightharpoonup press **CE** (the input is deleted).

Press Enter / Next



Weight

Fat %

Body Type

Weight
Healthy Fat %
Fat %
Body Type

Net
Weight Lock
Weight Lock Mode
Stabilized

Fat
Weight Lock Mode
Stabilized

Male/Female

The "step on" lamp flashes.

Note

• When **CE** is pressed, it returns to the previous screen.

Weight

Healthy Fat %

Fat %

Body Type

Standard | Athletic

24

How to use (Operating Instructions)

Step on the electrodes with bare feet.


Take off your socks and stockings before stepping on.



When the weight becomes stable, it changes to the screen at the right.

Note

- Do not step off the platform.
- This screen is not displayed if OFF is set in the "setting with or without an ID" (page 15). (The "body type selection" screen is displayed.)



Inpu


Input the ID number.

Input it by pressing 0 - 9

Note

- This screen is not displayed if OFF is set in the "setting with or without an ID" (* page 15).
- The ID number can be input from 0 9999999999. If Enter / Next is pressed, the non-inputted digits are filled with 0s.
- If it is mistakenly input,


 ⇒press CE (the input is deleted).
- If <u>CE</u> is pressed in the state that an ID number is not input, it returns to the "measurement start" screen.



Press Enter / Next

Note

• This screen is not displayed if OFF is set in the "setting with or without an ID" (page 15).



The lamp flashes on the "Body type"



Select the body type.

Press the body type selection keys to input. When the body type is selected, the lamp flashes on the "Gender".

Note

- This screen is not displayed if OFF is set in the "setting ON or OFF of the athletic mode selection" (Fig. page 16).
- The body type can also be selected with the numeric keys (1 2).
- If it is mistakenly input,
- ⇒press **CE** (the input is deleted, and it returns to the "body type selection" screen).
- If **CE** is pressed in the state that the body type is not input, it returns to the "ID number input" screen (or "measurement start" screen).

Weight

Select gender.

Press the male / female selection keys to input. When male or female is selected, the lamp flashes on the "Age".

Note

- If it is mistakenly input,
- ⇒press CE (the input is deleted, and it returns to the "gender selection" screen).
- If **CE** is pressed in the state that the male or female is not selected, it returns to the "body type selection" screen.



26

Male/Female

Input the age.

Input it by pressing 0 - 9

Note

- The age can be input from 5 99.
- If it is mistakenly input, \Rightarrow press **CE** (the input is deleted).
- If **CE** is pressed in the state that the age is not input, it returns to the "gender selection" screen.

Weight Male/Female Body Type

Press Enter / Next).


When the age is input, the lamp flashes on the




Input the height.

Input it by pressing 0 – 9

- The height can be input from 90.0 249.9 (90 249).
- If it is mistakenly input,
- ⇒press **CE** (the input is deleted).
- If **CE** is pressed in the state that the height is not input, it returns to the "age selection" screen.



Press Enter / Next



Under measurement of the body composition.

RAPPA display goes off sequentially.



Measurement completion

The measurement result and the body fat percentage evaluation are displayed.

It is automatically printed out. (In the case that other than 0 is set in the "Setting the number of sheets to print" (pp page 13).

Note


• Fat % Healthy range is not displayed if OFF is set in the "setting ON or OFF of the Fat % Healthy range display" (page 15).


Step off the platform.

It returns to the "measurement start" screen.



How to use (Operating Instructions)





How to use (Operating Instructions)

In the case of the one step mode

In the one step mode, after inputting the personal data, the weight and body composition are measured. **Body composition measurement** Personal data input Weight measurement

Press o/o to turn on the power.



Weight

∭ kg

Check that the body composition monitor is selected.

input the clothes weight.

Input it by pressing 0 – 9.


- The clothes weight (preset tare) can be input in the Body Type range of 0.0 - 10.0 kg.
- If it is mistakenly input, \Rightarrow press (CE) (the input is deleted).

Press (Enter / Next)



Fat % Body Type

• This screen is not displayed if OFF is set in the "setting with or without an ID" (page 15). (The "body type selection" screen is displayed.)



Input the ID number.

Input it by pressing () - (9)

Note

- This screen is not displayed if OFF is set in the "setting with or without an ID" (regrapage 15).
- The ID number can be input from 0 9999999999. If **Enter / Next** is pressed, the non-inputted digits are filled with 0s.
- If it is mistakenly input,
- \Rightarrow press **CE** (the input is deleted).
- If CE is pressed in the state that an ID number is not input, it returns to the "measurement start" screen.

Select the body type.

Press the body type selection keys to input. When the body type is selected, the lamp flashes on the "Gender".


- This screen is not displayed if OFF is set in the "setting ON or OFF of the athletic mode selection" (page 16).
- The body type can also be selected with the numeric keys (1 2)
- If it is mistakenly input,
- press (CE) (the input is deleted, and it returns to the "body type selection" screen).
- If CE is pressed in the state that the body type is not input, it returns to the "ID number input" screen (or "measurement start" screen).


Select gender.

Press the male / female selection keys to input. When male or female is selected, the lamp flashes on the "Age".

Note

- If it is mistakenly input,
- press (CE) (the input is deleted, and it returns to the "gender selection" screen).
- If CE is pressed in the state that the male or female is not selected, it returns to the "body type selection" screen.





How to use (Operating Instructions)

7 Ir

Input the age.

Input it by pressing **0** – **9**.

Note

- The age can be input from 5 99.
- If it is mistakenly input,
 ⇒ press CE (the input is deleted).
- If **CE** is pressed in the state that the age is not input, it returns to the "gender selection" screen.

Weight Lock Weight Lock Mode Stabilized


Fat %

Body Type

Standard Athletic

Press Enter / Next

When the age is input, the lamp flashes on the "Height".




Input the height.

Input it by pressing 0 - 9.

Note

- The height can be input from 90.0 249.9 (90 249).
- If it is mistakenly input,
- \Rightarrow press **CE** (the input is deleted).
- If **CE** is pressed in the state that the height is not input, it returns to the "age selection" screen.




Press Enter / Next



The "step on" lamp flashes.

Note

• When **CE** is pressed, it returns to the previous screen.



Step on the electrodes with bare feet.

Take off your socks and stockings before stepping on.



Under measurement of the body composition.

BBBB display goes off sequentially.



Measurement completion

The measurement result and the body fat percentage evaluation are displayed.

It is automatically printed out. (In the case that other than 0 is set in the "Setting the number of sheets to print" (page 13).

Note

• Fat % Healthy range is not displayed if OFF is set in the "setting ON or OFF of the Fat % Healthy range display" (page 15).


Step off the platform.


It returns to the "ID input" screen.

Note

• The "body type selection" screen is displayed in the case that OFF is set in the "setting with or without an ID" (page 15).







(GB)

How to use (Operating Instructions)

Input the target body fat ratio.

Input it by pressing 0 - 9

Note

- This screen is displayed if ON is set in the "setting ON or OFF of the target body fat ratio".
- The target body fat ratio can be input from 4-55.
- If it is mistakenly input,

press (CE) (the input is deleted).



If the number of sheet to print is set to "0", target body fat ratio setting function will be OFF, automatically.

If the target body fat is set to 0 or nothing, the target body fat ratio will not print.



Before you start a body weight management program and set the appropriate personal Caution body fat ratio, please consult your doctor. Tanita is not responsible for setting the appropriate target body fat ratio for specific individuals.


Press **O/O** to turn on the power.




Check that the body composition monitor is selected and input the clothes weight.

Input it by pressing 0 - 9 and .


- The clothing weight (tare) can be input in the range of 0.0 – 10.0 kg.
- If it is mistakenly input, \Rightarrow press **CE** (the input is deleted).



Press (Enter / Next)



• This screen is not displayed if OFF is set in the "setting with or without an ID" (page 15). (The "body type selection" screen is displayed.)



How to use (Operating Instructions)

Input an ID number.

Input it by pressing 0 - 9

Note

- This screen is not displayed in the case that OFF is set in the "setting with or without an ID" (page 15).
- The ID number can be input from 0 9999999999.
- \bullet If it is mistakenly input,
- ⇒press **CE** (the input is deleted).
- If CE is pressed in the state that an ID number is not input, it returns to the "measurement start."

Body Type

5

Press Enter / Next

Note


• This screen is not displayed in the case that OFF is set in the "setting with or without an ID" (page 15).



ID ► Net

Weight Lock Mode


Male/Female



The lamp flashes for "step on."

Note

• If **CE** is pressed, it returns to the previous screen.



6


Measurement completion

It is automatically printed out. (In the case that other than 0 is set in the "Setting the number of sheets to print" (page 13).



Step off the platform.

It returns to the "measurement start" screen.



Various criteria

- Criteria based on body fat percentage

Body fat percentage is the amount of body fat as a proportion of your body weight.

Reducing excess levels of body fat has shown to reduce the risk of certain conditions such as high blood pressure, heart disease, diabetes and cancer. The chart below shows the healthy ranges for body fat.

Body Fat Ranges for Standard Children 1 Body Fat Ranges for Standard Adults 2

How to use (Various criteria)

"New Body Fat Reference Curves for children"

¹ Susan Jebb et al. Obesity Reseach[o2] 2004;12:A156-157 ² Gallagher D et al. Am J Clin Nutr 2000,72:694-701 "Healthy percentage body fat ranges:an approach for developing guidelines based on body mass index."

		Und	ertat												Heal	lthy							0ve	ertat			Obe	ise																	
	5	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
	6	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
	7	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
	8	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
	9	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
	10	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
	11	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
	12	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
Female	13	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
Age	14	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
	15	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
	16	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
	17	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
	18	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
	19	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
	20-39	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	43	44	45
	40-59	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
	60-	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38 3	9 4	0 41	42	2 43	44	45
		٠,								-	00/										00/									_	00/									400	,				

	0	%								1	0%									2	20%									31	0%									40)%				
	5	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	6	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	7	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	8	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	9	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	10	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	11	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	12	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
Male	13	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
Age	14	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	15	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	16	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	17	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	18	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	19	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	20-39	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	40-59	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
	60-	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44 45
		Und	derfat											Hea	lthv											Ove	rfat				Ohes	e													

Copyright (C) 2004 TANITA Corporation, All Rights Reserved

Healthy Range Indicator

Your Body Composition Monitor automatically compares your body fat percentage reading to the Healthy Body Fat Range chart. After your body fat percentage has been calculated, a black bar will flash along the bottom of the display, identifying where you fall within the Body Fat Ranges for your age and gender.



- (-) : Underfat; below the healthy body fat range. Increased risk for health problems.
- (0): Healthy; within the healthy body fat percentage range for your
- (+) : Overfat; above the healthy range. Increased risk for health problems.
- (++): Obese: high above the healthy body fat range. Greatly increased risk of obesity-related health problems.

*Note: If you select Athlete mode, the unit will not display the Healthy Range Indicator.

> Athletes may have a lower body fat range depending on their particular sport or activity.

- What is total body water percentage?

Total Body Water Percentage is the total amount of fluid in a person's body expressed as a percentage of their total weight. Water plays a vital role in many of the body's processes and is found in every cell, tissue and organ. Maintaining a healthy total body water percentage will ensure the body functions efficiently and will reduce the risk of developing associated health problems.

Your body water levels naturally fluctuate throughout the day and night. Your body tends to be dehydrated after a longnight and there are differences in fluid distribution between day and night. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may cause variations in your hydration levels.

Your body water percentage reading should act as a guide and should not be used to specifically determine your absolute recommended total body water percentage. It is important to look for long-term changes in total body water percentage and maintain a consistent, healthy total body water percentage.

Drinking a large quantity of water in one sitting will not instantly change your water level. In fact, it will increase your body fat reading due to the additional weight gain. Please monitor all readings over time to track the relative change.

Every individual varies but as a guide the average total body water percentage ranges for a healthy adult are:

Female: 45 to 60% **Male** : 50 to 65%

Source: Based on Tanita's Internal Research

Note: The total body water percentage will tend to decrease as the percentage of body fat increases. A person with a high percentage of body fat may fall below the average body water percentage. As you lose body, fat the total body water percentage should gradually move towards the typical range given

- What is visceral fat rating?

This feature indicates the rating of visceral fat in your body.

Visceral fat is the fat that is in the internal abdominal cavity, surrounding the vital organs in the trunk (abdominal) area. Research shows that even if your weight and body fat remains constant, as you get older the distribution of fat changes and is more likely to shift to the trunk area especially post menopause. Ensuring you, have healthy levels of visceral fat my reduce the risk of certain diseases such as heart disease, high blood pressure, and the onset of type 2 diabetes.

The Tanita Body Composition Monitor will provide you with a visceral fat rating from 1 – 59.

Rating from 1 to 12

Indicates you have a healthy level of visceral fat. Continue monitoring your rating to ensure that it stays within the healthy range.

Rating from 13 – 59

Indicates you have an excess level of visceral fat. Consider making changes in your lifestyle possibly through diet changes or increasing exercise.

Source: Data from Columbia University (New York) & Tanita Institute (Tokyo)

- Even if you have a low body fat rate, you may have a high visceral fat level.
- For medical diagnosis, consult a physician.

- What is basal metabolic rate (BMR)?

WHAT IS BMR?

Your Basal Metabolic Rate(BMR) is the minimum level of energy your body needs when at rest to function effectively including your respiratory and circulatory organs, neural system, liver, kidneys, and other organs. You burn calories when sleeping.

About 70% of calories consumed every day are used for your basal metabolism. In addition, energy is used when doing any kind of activity however; the more vigorous the activity is the more calories are burned. This is because skeletal muscle (which accounts for approximately 40% of your body weight) acts as your metabolic engine and uses a large amount of energy. Your basal metabolism is greatly affected by the quantity of muscles you have, therefor increasing your muscle mass will help increase your basal metabolism.

By studying healthy individuals, scientists have found that as people age, their metabolic rate changes. Basal metabolism rises as a child matures. After a peak at the age of 16 or 17, it typically starts to decrease gradually.

Having a higher basal metabolism will increase the number of calories used and help to decrease the amount of body fat. A low basal metabolic rate will make it harder to lose body fat and overall weight.

HOW DOES A TANITA BODY COMPOSITION MONITOR CALCULATE BMR?

The basic way of calculating Basal Metabolic Rate BMR is a standard equation using weight and age. Tanita has conducted in-depth research into the relationship of BMR and body composition giving a much more accurate and personalized reading for the user based on the impedance measurement. This method has been medically validated using indirect calorimetry (measuring the breath composition).*

* Reliability on equation for Basal Metabolic Rate: At 2002 Nutrition Week : A Scientific and Clinical Forum and Exposition Title: International Comparison: Resting Energy Expenditure Prediction Models: The American Journal of Clinical Nutrition

- What is metabolic age?

This feature calculates your BMR and indicates the average age associated with that type of metabolism. If your BMR Age is higher than your actual age, it is an indication that you need to improve your metabolic rate. Increased exercise will build healthy muscle tissue, which will improve your metabolic age.

You will obtain a reading between 12 and 50. Under 12 will be displayed as "12" and over 50 displayed as "50".

- What is mascle mass?


This feature indicates the weight of muscle in your body, the muscle mass displayed includes the skeletal muscles, smooth muscles (such as cardiac and digestive muscles) and the water contained in these muscles.

Muscles play an important role as they act as an engine in consuming energy. As your muscle mass increase, your energy consumption increases helping you reduce excess body fat levels and lose weight in a healthy way.

- What is physique rating?

This feature assesses your physique according to the ratio of body fat ad muscle mass in your body. As you become more active and reduce the amount of body fat, your physique rating will also change accordingly. Even though your weight may not change, your muscle mass and body fat levels may be changing making you healthier and at lower risk of certain diseases.

Each person should set their own goal of which physique they would like and follow a diet and fitness program to meet that goal.




(GB)

How to use (Various criteria)

(GB)

- Please check the following before asking for repair.

Please check **Symptom Impedance** measurement error • Measure with bare feet. ErryO • When the soles of your feet are dry, drop water for about 0.5 ml with the attached dropper on the electrodes before measurement. • Check the input information. Err67 How to **Zero point error** measure • Turn off the power, and remove the items on the platform and turn on the power again, and then redo the measurement. • Is it installed at a place with vibrations? • Is the platform inclined? ⇒Keep the platform horizontal. (☐ pages 6 The measured and 7) weight is not stable. • Is anything caught in the gaps of the platform? Remove anything caught in the gaps. Nothing is displayed even after • Check that the power supply is connected turning on the correctly. Display power. part is • The weight to measure exceeds the measurement range. displayed.



• In the case that an error is displayed other than the above, turn off the power once, and then measure again.

If the same error is displayed repeatedly, contact our customer service centre.

Connection with a personal computer

B)

When necessary (Connection with a personal computer)




The RS-232C interface enables input and output from this equipment. This equipment is not capable of being remotely controlled by external equipment, such as a computer.

- Specifications

Communication standards	Compatible with EIA RS-232C
Communication method	Asynchronous communication method
Signal speed	9600 bps
Data bit length	8 bits
Parity	NONE
Stop bit	1 bit
Flow control	NONE
Terminator	CR+LF

- Signal Names and Connection Methods


Terminal Number	Signal Name
1	※ 1
2	RXD
3	TXD
4	※ 1
5	GND
6	※ 1
7	※ 2
8	% 2
9	No Connection



*1: Pin Nos. 1, 4 and 6 are internal connections.

- Connection Example

Please be sure to use a straight cable when the equipment is connected to an external computer.





A modem cable cannot be used.

^{*2:} Pin Nos. 7 and 8 are internal connections.

When necessary (Connection with a personal computer)

Transmission Data is output immediately after measurement regardless of the status of Transmission Data is output immediately after measurement regardless of the status of the receiving equipment (personal computer, etc.). Therefore the receiving equipment needs to be ready to accept the data before measuring.

- Transmit data

Transmit data is output immediately after measurement regardless of the state of the receiving side (personal computer, etc.). Therefore, the receiving side must always be ready to receive before measurement.

*PC mode is a mode to send personal data from the personal computer side and to receive measurement results.

(1) Output data format

The measured data is output in the following format.

- Each piece data is with a comma-delimited (,).
- The terminator (the end of the data) is CR (ASCII code 0DH), LF (ASCII code 0AH). 0 herein is zero.

Entire body data

N	/lodel	Se	rial No.	I D	number	Date (d	d/mm/yyyy)	T	ime
MO	"XXXXXX"	SN	"XXXXXXXX"	ID	"XXXXXXXXXX"	Da	"dd/mm/yyyy"	TI	"hh:mm"
Boo	dy type	G	ender		Age	H	leight	Cloth	es (tare)
Bt	Oor2	GE	1or2	AG	XX	Hm	XXX.X	Pt	XX.XX
1/1	/eight	Por	dy fat %	E _r	at mass	Eat f	ree mass	Mued	ele mass
Wk	XXX.X	FW	XX.X	fW	XXX.X	MW	XXX.X	mW	XXX.X
									<u>'</u>
Mus	cle score	Bor	ne mass		TBW	T	3W %	E	3MI
sW	XX	bW	XXX.X	wW	XXX.X	WW	XXX.X	MI	XXX.X
	1				1				
Standard	body weight	Degree	of obesity	Viscer	al fat rating	BN	/IR (kJ)	BMF	R (kcal)
Sw	XXX.X	OV	XX.XX	IF	XX	rb	XXXXX	rB	XXXXX
	_		1						
BM	R score	Meta	bolic age	Rohr	er's index	Target	body fat %	Predict	ed weight
rJ	XX	rA	XX	R0	XXXX.X	gF	XX	gW	XXX.X
	i								
Predict	ed fat mass	Fat to	gain / lese	Im	pedance	Ch	ecksum		
gf	XXX.X	gt	XXXX.X	ZF	XXXX.X	CS	XX		

^{*} Each value is comma-delimited (,).

(2) Output data items

					Order o	f output	
Item	Header	Format	Contents	Body c	omposition	monitor	6
				Adult	Athlete	Child	Scale
Control data	{0	Fix to 16	2 byte fixed length	1	1	1	1
Control data	~0	Fix to 1	1 byte fixed length	2	2	2	2
Control data	~1	Fix to 1	1 byte fixed length	3	3	3	
Control data	~2	Fix to 1	1 byte fixed length	4	4	4	
Model	MO	"XXXXXX"	8 byte fixed length	5	5	5	3
Serial No.	SN	"XXXXXXXX"	10 byte fixed length	6	6	6	4
ID number	ID	"XXXXXXXXXX"	12 byte fixed length	7	7	7	5
Date (dd/mm/yyyy)	Da	"dd/mm/yyyy"	12 byte fixed length	8	8	8	6
Time	TI	"hh:mm"	7 byte fixed length	9	9	9	7
Body type	Bt	0 or 2	1 byte fixed length (0: standard 2: athlete)	10	10	10	
Gender	GE	1 or 2	1 byte fixed length (1: male 2: female)	11	11	11	
Age	AG	XX	"1 – 2 byte variable length (unit: age, right-aligned)"	12	12	12	
Height	Hm	XXX.X	"4 – 5 byte variable length, to 1 place of decimals (unit: cm)"	13	13	13	
Clothes (tare)	Pt	XX.X	"3 – 4 byte variable length, to 1 place of decimals (unit: kg)"	14	14	14	8
Weight	Wk	XXX.X	"3 – 5 byte variable length, to 1 place of decimals (unit: kg)"	15	15	15	9
Body fat %	FW	XX.X	"3 – 4 byte variable length, to 1 place of decimals (unit: %)"	16	16	16	
Fat mass	fW	XXX.X	"3 – 5 byte variable length, to 1 place of decimals (unit: kg)"	17	17	17	
Fat free mass	MW	XXX.X	"3 – 5 byte variable length, to 1 place of decimals (unit: kg)"	18	18	18	
Muscle mass	mW	XXX.X	"3 – 5 byte variable length, to 1 place of decimals (unit: kg)"	19	19	19	
Muscle score	sW	XX	1 – 2 byte variable length (1-24)	20	20		
Bone mass	bW	XXX.X	"3 – 5 byte variable length, to 1 place of decimals (unit: kg)"	21	21		
TBW	wW	XXX.X	"3 – 5 byte variable length, to 1 place of decimals (unit: kg)"	22	22	20	
TBW %	ww	XXX.X	"3 – 5 byte variable length, to 1 place of decimals"	23	23	21	
BMI	MI	XXX.X	"3 – 5 byte variable length, to 1 place of decimals (unit: kg)"	24	24		
Standard body weight	Sw	XXX.X	"3 – 5 byte variable length, to 1 place of decimals (unit: kg)"	25			
Degree of obesity	OV	XX.XX	"3 – 5 byte variable length, to 1 place of decimals (unit: %)"	26			
Visceral fat rating	IF	XX	1 – 2 byte variable length	27	25		
BMR (kJ)	rb	XXXXX	1 – 5 byte variable length (unit: kJ)	28	26		
BMR (kcal)	rB	XXXXX	1 – 5 byte variable length (unit: kcal)	29	27		
BMR score	rJ	XX	1 – 2 byte variable length	30	28		
Metabolic age	rA	XX	2 byte fixed length	31	29		
Rohrer's index	RO	XXXX.X	4 – 6 byte variable length			22	
Target body fat %	gF	XX	1 – 2 byte variable length	32	30	23	
Predicted weight	gW	XXX.X	3 – 5 byte variable length	33	31	24	
Predicted fat mass	gf	XXX.X	3 – 5 byte variable length	34	32	25	
Fat to gain / lese	gt	XXXX.X	3 – 6 byte variable length	35	33	26	
Impedance	ZF	XXXX.X	5 – 6 byte variable length	36	34	27	
Checksum	CS	XX	2 byte fixed length	37	35	28	10

Technical notes

(GB)

When necessary (Technical notes)

Introduction

Body composition measurement by the BIA method.

This equipment provides estimated values for each measured value of body fat percentage, fat mass, fatfree mass, muscle mass and bone mass by the DXA method, estimated value for the total body water measured value by the dilution method and estimated value for the visceral fat rating by MRI method using the Bioelectrical Impedance Analysis (BIA method).

For measurement, a mode must be selected based on body type.

- 1) Standard (for 5-99 years of age)
- 2) Athletic (for Athletic persons who exercise considerably more than non-athlete)

Making a distinction by body type in the measurement mode produces more reliable body composition measurements for athletic persons, whose body compositions differ from those of average persons.


- Principles of body composition measurement

BIA is a means of measuring body composition - fat mass, predicted muscle mass, etc. - by measuring bioelectrical impedance in the body. Fat within the body allows almost no electricity to pass through, while electricity passes rather easily through water, much of which is found in muscles. The degree of difficulty with which electricity passes through a substance is known as the electrical resistance, and the percentage of fat and other body constituents can be inferred from measurements of this resistance.

The Tanita Body Composition Analyzer measures body composition using a constant current source with a high frequency current (50kHz, 90μA). The 8 electrodes are positioned so that electric current is supplied from the electrodes on the tips of the toes of both feet, and voltage is measured on the heel of both feet. The current flows into the upper limbs or lower limbs, depending on the body part(s) to be measured.

- What is the DXA method?

DXA was originally designed to measure bone mineral content, but in the full-body scan mode the body fat percentage, fat mass, and fat free mass of individual body parts (arms, legs, trunk) can also be measured. The image below shows one example of body composition measurement results obtained by DXA.



obtained by DXA (Lunar Co., Ltd; DPX-L)

- What is dilution method?

In the dilution method, a labeled substance for a known amount is given and the concentration in equilibrium diffusing evenly is measured to obtain the total amount of the solvent that dilutes the labeled substance.

To measure the total body water (TBW), deuterium oxide (D2O) is generally used as the labeled substance. Deuterium oxide uses the overall total body water as dilution space so the total body water can be obtained. To obtain the extracellular fluid amount, sodium bromide (NaBr) is used as a labeled substance. Bromine (Br) is said to not enter the inside of cells, and uses extracellular fluid as the dilution space.


- What is the visceral fat?


Visceral adipose tissue (VAT) is fat that accumulates in the abdominal cavity and around internal organs. VAT is said to be more likely to cause lifestyle-related diseases than subcutaneous adipose tissue (SCAT). Accordingly, knowing and periodically checking the VAT accumulation risk serves as an important guide in the prevention of lifestyle-related diseases.


Tanita has developed the technology for measuring the VAT accumulation risk through bioelectrical impedance analysis (BIA) in comparison with image analysis applied to magnetic resonance imaging (MRI), in addition to the established technology for measuring the percent of body fat. The VAT accumulation risk is calculated by estimating the VAT area by the BIA method on the basis of MRI image processing. This method has a higher correlation than the estimation of the VAT accumulation risk based on BMI or abdominal circumference (waist circumference), allowing estimation that corresponds more precisely to individuals.

*The VAT area by MRI is calculated by carrying out an image processing of the cross section of the lumber vertebra L4-L5 regions.


(Fig. 1 - Fig. 3: Research results by N. Y. Columbia University and Jikei University Published by the North American Association for the Study of Obesity [NAASO] in 2004.)








<Fig. 3> Relationship between VAT Area by MRI and Estimated VAT Area by Tanita's BIA





In the BIA method, impedance is measured and the body composition is calculated based on the value. It is known that impedance changes by the amount of the total body water that occupies about 60% of weight and the change in its distribution and temperature change. Therefore, for the purpose of research or for daily repeating of measurements, the measurement conditions must be kept constant. Measurement under the changing conditions of temperature and total body water distribution or blood flow volume of extremities due to exercising, taking a bath, etc., affects the measurement result since the electric resistance in the body also changes.

Therefore, it is recommended to measure under the following conditions for stable measurement.

- 1) 3 hours have passed after getting up and normal lifestyle activities are carried out during this period. (The impedance transits staying at a high level if you remain sitting after getting up or drive a car, etc.)
- 2) 3 hours or more have passed after eating. (For 2 3 hours after eating, the impedance has a tendency to decrease.)
- 3) 12 hours or more have passed after vigorous exercise for measurement. (The tendency toward changes in impedance is not stable depending on the type and rigorousness of the exercise.)
- 4) If possible urinate before taking measurement.
- 5) For repeated measurements, measure at the same hour as much as possible. (At the same time of measurement of weight, the measurements can be made more stable by measuring at the same time of the day)

Very stable measured values can be obtained by measuring under the above conditions.

And in the development of this equipment, the following 6 items were set as conditions for the regression equation.

- 1) Prohibition of alcohol intake for 12 hours before measurement
- 2) Prohibition of excessive exercise for 12 hours before measurement.
- 3) Prohibition of excessive eating and drinking the day before measurement
- 4) Prohibition of eating and drinking for 3 hours before measurement
- 5) Avoidance of the menstrual period (women)


2) Inter-day changes


The diagrams below offer examples of actual measurements made of inter-day changes. A study was done to determine the degree of change in the impedance between the feet during dehydration; the first two days represent a normal daily routine, while in the latter two days a state of dehydration was induced using a sauna.

No significant inter-day change was measured in body weight, impedance between the feet, or body fat percentage during the normal daily routine. During the dehydrated state, however, a drop in body weight of 1kg was noted, with the impedance between the feet rising approximately $15\,\Omega$ on the first day of dehydration and $30\text{-}35\,\Omega$ on the second day. As a result, body fat percentage was up by around 1% on the first day of dehydration and by 1.5-2% on the second day.

As mentioned earlier, impedance increases when body weight is reduced (such as by dehydration), and decreases when body weight is increased through excess consumption of food and drink. The inter-day change in impedance is thus inversely proportional to the change in body weight.

These inter-day changes stem from such causes as:





- 1) Temporary increases in body weight (total body water) through overeating and overdrinking
- 2) Dehydration due to heavy sweating during vigorous exercise
- 3) Dehydration due to alcohol consumption or the use of diuretics
- 4) Dehydration due to heavy sweating during saunas, etc.

Accordingly, it is recommended that instructions be provided to the subject to help eliminate these causes when accurate measurements are needed.

(GB)

When necessary (Technical notes)


The New Regression Formula for Basal Metabolic Rate (BMR)


It has long been said among medical and nutritional specialists that "The Basal Metabolic Rate (BMR) is more determined by the Fat Free Mass (FFM) than by the body weight" (Persons of a given body weight with a higher FFM will have a higher BMR), and that from the aspect of evaluating the body composition, should be estimated from the FFM. In addition, in cases of simple estimation formulae which can calculate from the height, weight and age, without evaluating the body composition, there was a problem with excessively high evaluations being given to obese persons with large body weight, and conversely excessively small BMR evaluations given to muscular athletes, though these are not as many in number. Currently, the BMR estimation recursion formula developed by Tanita, the manufacturer of body composition analyzers, based on their research, works by multiple regressive analysis using this FFM, and has a higher degree of accuracy in the individual differences in body composition. In order to derive the BMR, resting respiratory metabolism (Resting Energy Expenditure: REE) was measured using a breath gas analysis device, and this estimation recursion formula was created based on this data.

<Figure 1> The Relationship Between Resting Energy Expenditure (REE) According to Breath Gas Analysis and Weight, FFM

(Presented at Nutrition Week, Held in San Diego in 2002)

As shown in Figure 1: the REE (BMR) has a stronger relationship to the FFM than to body weight, and a difference is visible between males and females in the distribution trends. We see that in principle that we should calculate from the FFM rather than by the old formula centered on the relationship with weight.






<Figure 2> Comparison of BMR Values from the TANITA Multiple regression model and Breath Analysis

(Presented at Nutrition Week Held in San Diego in 2002)

The current BMR retrogression formula is a formula which acts on the principle of using the FFM value from the results of body composition measurement according to the BIA, A good relationship is shown in the BMR value based on actual breath analysis REE or R=0.9 (p<0.0001). These results were presented at the First Annual Nutrition Week (American College of Nutrition, American Society for Clinical Nutrition, American Society for Parenteral and Enteral Nutrition, North American Association for the Study of Obesity) held in 2002 in San Diego.

NOTE: This model has been calibrated for those between ages of 18-84. Those individuals outside of this age range may not be obtain accurate readings.



Specifications

Model		SC-330
		AC adapter (included) Centre Minus MODEL SA165A-0950U-3 CLASS 2
Power source		Input Voltage: 100-240 VAC 50/60 Hz 1.5A Output Voltage: 7 VDC
		Rated Current: 4 A No Load Input Voltage: 7VDC
Power Consum	nption	28 W
	Measurement System	Tetra polar Bioelectrical Impedance Analysis
	Measurement Frequency	50 kHz
Immedance	Measurement Current	90 μΑ
Impedance Measurement	Electrode Material	Pressure Contact Stainless Steel Foot Pads
11204542 01110110	Measurement Style	Between Both Feet
	Measurement Range	150 - 1200Ω
	Accuracy at first calibration	± 2%
Weight	Measurement System	Strain Gauge Load Cell
Measurement	Maximum Capacity / Minimum Graduation	270 kg / 0.1 kg
	Accuracy at first calibration	± 0.2 kg
	Clothes Weight	0 - 10 kg / 0.1 kg increments
	Gender	Male / Female
Input Items	Body Type	Standard (5 - 99 years) / Athletic (18 - 99 years)
input items	Age	5 - 99 years old / 1 year increments
	Height	90 - 249.9 cm / 0.1 cm increments
	Target Body Fat %	4 - 55 %

Target Body Fat % Weight Gender Body Type Age Height FAT% Logo Model Name Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass BMR	4 - 55 % 0 - 270 kg / 0.1 kg increments Male / Female Standard / Athletic 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments TANITA LOGO (240x64 dot) SC-330 2005 / 1 / 1 - 2099 / 12 / 31 00000000 - 99999999 0000000000 - 99999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments 0.1 kg increments
Gender Body Type Age Height FAT% Logo Model Name Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	Male / Female Standard / Athletic 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 3 - 75% / 0.1% increments TANITA LOGO (240x64 dot) SC-330 2005 / 1 / 1 - 2099 / 12 / 31 00000000 - 999999999 0000000000 - 999999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 0 - 270 kg / 0.1% increments 0.1 kg increments
Body Type Age Height FAT% Logo Model Name Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	Standard / Athletic 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 3 - 75% / 0.1% increments TANITA LOGO (240x64 dot) SC-330 2005 / 1 / 1 - 2099 / 12 / 31 000000000 - 999999999 0000000000 - 999999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 0 - 270 kg / 0.1% increments 0.1 kg increments
Age Height FAT% Logo Model Name Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 3 - 75% / 0.1% increments TANITA LOGO (240x64 dot) SC-330 2005 / 1 / 1 - 2099 / 12 / 31 00000000 - 99999999 0000000000 - 999999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 0 - 10 kg increments 0.1 kg increments
Height FAT% Logo Model Name Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	90 - 249.9 cm / 0.1 cm increments 3 - 75% / 0.1% increments TANITA LOGO (240x64 dot) SC-330 2005 / 1 / 1 - 2099 / 12 / 31 00000000 - 999999999 0000000000 - 999999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 0 - 270 kg / 0.1% increments 0.1 kg increments
FAT% Logo Model Name Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	3 - 75% / 0.1% increments TANITA LOGO (240x64 dot) SC-330 2005 / 1 / 1 - 2099 / 12 / 31 00000000 - 99999999 0000000000 - 999999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 0 - 270 kg / 0.1% increments 0.1 kg increments 15 - 85% / 0.1% increments
Logo Model Name Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	TANITA LOGO (240x64 dot) SC-330 2005 / 1 / 1 - 2099 / 12 / 31 00000000 - 999999999 0000000000 - 999999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 3 - 75% / 0.1% increments 0.1 kg increments 15 - 85% / 0.1% increments
Model Name Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	SC-330 2005 / 1 / 1 - 2099 / 12 / 31 000000000 - 999999999 0000000000 - 9999999999
Date and Time Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	2005 / 1 / 1 - 2099 / 12 / 31 00000000 - 99999999 0000000000 - 999999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 0 - 270 kg / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Serial No. ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	00000000 - 999999999 000000000 - 9999999999
ID Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	000000000 - 9999999999 Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 3 - 75% / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Body Type Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	Standard (5 - 99 years) / Athletic (18 - 99 years) Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 3 - 75% / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Gender Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	Male / Female 5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 3 - 75% / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Age Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	5 - 99 years old / 1 year increments 90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 3 - 75% / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Height Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	90 - 249.9 cm / 0.1 cm increments 0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 3 - 75% / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Clothes Weight Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	0 - 10 kg / 0.1 kg increments 0 - 270 kg / 0.1 kg increments 3 - 75% / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Weight FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	0 - 270 kg / 0.1 kg increments 3 - 75% / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
FAT% Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	3 - 75% / 0.1% increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Fat Mass FFM Muscle Mass TBW TBW% Predicted Bone Mass	0.1 kg increments 0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
FFM Muscle Mass TBW TBW% Predicted Bone Mass	0.1 kg increments 0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
Muscle Mass TBW TBW% Predicted Bone Mass	0.1 kg increments 0.1 kg increments 15 - 85% / 0.1% increments
TBW TBW% Predicted Bone Mass	0.1 kg increments 15 - 85% / 0.1% increments
TBW% Predicted Bone Mass	15 - 85% / 0.1% increments
Predicted Bone Mass	
	0.1 kg increments
	1 kJ increments / 1 kcal increments
Metabolic Age	1 years increments (12 - 90 years)
Visceral FAT Rating	1 Level increments (1-59 Level)
BMI	0.1 increments
Predicted Weight	0.1 kg increments
Desirable range	3 - 75% / 0.1% increments
FAT % Graph	5 - 75% / 0.1% increments
BMI Graph	
	150 - 1200Ω
Impedance	3 Rows, 5 Digits LCD
	RS-232C (D-sub 9 pins Female Connector)
100	0 - 35°C
150	30 - 80% (without condensation)
note Display Version)	6.8 kg
	12.1kg
	1w.1ng
Platform	372 × 375 × 101 mm
r	Visceral FAT Level Graph Muscle Mass Graph BMR Graph Physique Rating Impedance mote Display Version) lumn Mounted Version)

54

When necessary (Specifications)

(E

This device features radio interference suppression in compliance with EC Regulation 89/336/EC

TANITA EUROPE B.V. Holland Office Centre, Kruisweg 813-A 2132NG Hoofddorp, the Netherlands Tel: +31 (0) 23-5540188 FAX: +31 (0) 23-5579065 http://www.tanita.eu

TANITA UK LTD.

The Barn, Philpots Close, Yiewsley, Middlesex, UB7 7RY, United Kingdom Tel: +44 (0) 1895-438577 FAX: +44 (0) 1895-438511 http://www.tanita.co.uk

TANITA Corporation 1-14-2, Maeno-cho, Itabashi-ku, Tokyo, Japan Tel: +81 (0) 3-3968-2123 / +81 (0) 3-3968-7048 FAX: +81 (0) 3-3968-7048 FAX: +81 (0) 3-3968-7048 FAX: +81 (0) 3-3968-7048

TANITA Corporation of America, Inc.

2625 South Clearbrook Drive Arlington Heights, Illinois 60005, USA Tel: +1 847-640-9241 FAX: +1 847-640-9261 http://www.tanita.com

TANITA Health Equipment H.K.LTD.

Unit 301-303, 3/F Wing On Plaza, 62 Mody Road, Tsimshatsui East, Kowloon, Hong Kong Tel: +852 2838-7111 FAX: +852 2838-8667